


INSIDE CONFIGURATION



SPECIFICATIONS

ITEM	QUANTITY	DESCRIPTION
DO Down Tube	1	Down Tube for the Non-Invasive DO Sensor Port, Fits Standard 12mm Probes
Harvest Down Tube	1	Down Tube for the Harvest Line, 1/4" ID
Impeller	1	Pitched Blade Impeller, Polycarbonate
Microsparger Down Tube	1	Down Tube for the Microsparger, 1/8" ID
Optical pH Down Tube	1	Down Tube for the Non-Invasive Optical pH Port, 1/4" ID
Sample Down Tube	1	Down Tube for the Sample Port. 1/8" ID
Shaft	1	Impeller Shaft, Polycarbonate
Thermowell / RTD Down Tube	1	Down Tube for the Thermowell / RTD Port, Accommodates up to 1/4" RTD